

	Nursing Anne Simulator (Sim-Pad Capable) 11-81-1505	Nursing Anne Manikin (Sim-Pad Capable) 11-81-0505	Nursing Anne Manikin (Non-SimPad Capable) 11-81-0050
GENERAL			
Latex-Free	✓	—	—
Available in ethnic skin tones	✓ 23	✓ 22	✓ 22
BIRTH MONITORING			
Fetal heart tones	✓ 6	✓ 4	—
SIMSTORE COURSEWARE			
Courseware available at SimStore.com	✓	✓	—
EDUCATIONAL SERVICES (OPTIONAL)			
Introduction course	✓	✓	—
Customized introduction	✓	✓	—
Customized scenario programming	✓	✓	—
TECHNICAL SERVICES (OPTIONAL)			
Installation	✓	—	—
Extended warranty	✓	—	—
Preventative maintenance	✓	—	—
ValuePlus Package	✓	—	—
Loaner Program (Americas only)	✓	✓	—
AIRWAY FEATURES			
Intubation	✓	✓ 3	✓ 3
Oral/Nasal airway insertion	✓	✓	✓
Endotracheal tubes – insertion, securing, and care	✓	✓	✓
Oropharyngeal and nasopharyngeal airways – insertion and suctioning	✓	✓	✓
Supraglottic airways (LMA, LTA, and Combitube)	✓	—	—
Right mainstream intubation	✓	—	—
Sellick maneuver	✓	—	—
Oxygen delivery procedures	✓	✓	✓
Suctioning techniques	✓	✓	✓
Tracheostomy care and suctioning	✓	✓	✓
Trach opening for placement	✓	✓	✓
Auscultation of lung sounds	✓ 9	✓ 4	—
Auscultation of lung sounds during ventilation	✓ 6	—	—
Lung sounds, synchronized with breathing rate	✓ 6	✓ 4	—
Individual lung or bilateral sound selection	✓ 6	✓ 4	—
Laryngeal mask airway	✓	—	—
Positive pressure ventilation	✓	—	—
Realistic chest rise and fall	✓	—	—
CIRCULATION			
Instructor-controlled ECG rhythm	✓ 3	✓ 4	—
ECG monitoring	—	✓ 4	—
Auscultation of heart sounds	✓ 6	✓ 4	—
Heart sounds, synchronized with programmable ECG	✓ 6	✓ 4	—
BLOOD PRESSURE/PULSES			
Adjustable pulse strength	✓ 6	✓ 4	—
Adjustable heart rate	✓ 6	✓ 4	—
Automatic carotid pulse (electronic pulses)	✓ 7	—	—
Synchronized (with ECG)	✓ 6	✓	✓
Manual carotid pulses (pulse bulb)	—	✓	✓
Carotid pulse check sensors	✓ 6	—	—
Brachial pulses	✓ 7	✓	—
Pulses only active when palpated	✓ 6	✓ 4	—
Pulse strengths depended on blood pressure	✓ 6	✓ 4	—
Adjustable blood pressure levels (systolic/diastolic)	✓ 6	✓ 4	—
Blood pressure auscultation	✓ 7	✓ 4	—
Korotkoff sounds, synchronized with ECG	✓ 6	✓ 4	—
Auscultative gap, with on/off feature	✓ 6	✓	✓
Blood pressure/pulses	✓ 7	—	—
Femoral pulse (electronic)	✓ 7	—	—
Dorsalis pedis (pedal) pulse (electronic)	✓ 7	—	—
Radial pulse (electronic)	✓ 7	—	—
CATHETERIZATION			
Complete urinary catheterization	✓	✓	✓
DEBRIEFING			
SimView	✓ 17	✓ 17	—
Video debriefing	✓ 17	✓ 17	—
Debriefing through recorded events	✓ 6	✓ 4	—
SimView Mobile	✓ 17	✓ 17	—
BREATHING			
Spontaneous breathing without external compressor	✓	—	—
CIRCULATORY SKILLS AND IV DRUG ADMINISTRATION			
IM Injections	✓	✓	✓
Volume infusion	✓	✓	✓
Realistic flashback	—	✓	✓
Central catheters – insertion, securing, dressing, and maintenance	✓	—	—
Catheter openings – subclavian, jejunostomy, and Hickman (placement)	—	✓	✓
Subcutaneous heparin and insulin injection – abdominal	✓	—	—
IV insertion	✓ 21	✓ 20	✓ 20
Preported bilateral IV access	✓	—	—
Central Line – fluid infusion	✓	—	—
CPR			
Anatomical landmarks	✓	✓	✓
Ventilation with bag-valve-mask	✓	✓	✓
Chest compressions	✓	—	—
GASTROINTESTINAL			
Independent stomach	✓	✓	✓
NG Tubes – insertion, care, and removal	✓	✓	✓
Gastric lavage and gavage	✓	✓	✓
Gastronomy tube – placement and care	✓	—	—
Nasoenteric and esophageal tubes – insertion, care, and removal	✓	✓	✓
Feeding tubes – insertion and removal	✓	✓	✓
Auscultation of bowel sounds – normal/abnormal	✓ 6	✓ 4	—
Interchangeable stoma sites	—	✓ 10	✓ 10
Colostomy irrigation	✓	✓	✓
Rectal suppository administration	✓	—	—
Enema procedures with rectum	✓	✓	✓
Infected colostomy stoma	✓	—	—
Ostomy set: normal, non-perfusing (dusky) and infected	✓	—	—
Gastronomy tube – feeding	✓	—	—
EYE FEATURES			
Interchangeable pupils	✓	—	—
Unequal pupils	✓ 12	—	—
Dilated pupil	✓ 12	—	—
Blinking eyelids	✓	—	—
WOMEN'S HEALTH			
Vaginal douching	—	✓	✓
Breast examination	—	✓ 13	✓ 13
Recognition of breast disorders	—	✓ 13	✓ 13
Size of nodules and relative locations	—	✓ 13	✓ 13
Interchangeable abnormalities	—	✓ 13	✓ 13
Postoperative left total mastectomy	—	✓ 14	✓ 14
Fundus assessment and massage	—	✓ 15	✓ 15
Anal hemorrhoids	—	✓ 15	✓ 15
Episiotomy – 2nd degree midline	—	✓ 15	✓ 15
Perineal care	✓	—	—
Vaginal suppository/medication insertion	✓	—	—
GENDER/AGE/SIZE			
Adult	✓	✓	✓
Male	✓ 1	—	—
Female	✓	✓	✓
Full-body manikin	✓	✓	✓
EXTRICATION/IMMOBILIZATION			
Extrication	✓	—	—
Immobilization	✓	—	—
Victim handling	✓	—	—
SIMULATION ADMINISTRATION			
Event logging	✓ 6	✓ 4	—
Remote control	✓ 11	✓ 4	—
Create and edit scenarios	✓ 6	✓ 4	—
Pre-programmed scenarios	✓	✓ 4	—
Web camera	✓ 18	✓ 17	—
Simulated patient monitor	✓ 18	✓ 17	—
Preset patient state levels	✓ 6	✓ 4	—
Patient monitor with configurable vital sign display	✓ 18	✓ 17	—
PC control	✓	—	—
Instructor PC	✓ 17	—	—
Log	✓	✓ 19	—
Laptop PC	✓ 17	—	—
Wireless operation	✓	—	—
Create and edit themes	✓	✓ 19	—
VOCAL SOUNDS			
Pre-recorded vocal sounds	✓ 6	✓ 4	—
Live vocal sounds (wireless or wired mic)	✓ 8	✓ 4	—
WOUNDS/TRAUMA			
Nasal packing	✓	✓	✓
Optional wound care and assessment module set	✓	✓	✓
OTHER			
Realistic patient handling – full range of motion	✓	✓	✓
Eye/ear irrigation (simulated)	—	✓	✓
Application/instillation of medications – eye, ear, nose	—	✓	✓
Mouth care	✓	✓	✓
Denture care	✓ 2	✓	✓
SpO ₂ saturation	✓ 16	—	—
Sits unassisted	✓	—	—
Bends at waist	✓	—	—

- 1) Requires optional male genitalia
- 2) Upper
- 3) oral, nasal, digital
- 4) requires optional SimPad System
- 5) with SimPad Classic, SimPad PLUS or LLEAP and SimPad/LLEAP Patient Monitor
- 6) with SimPad Classic, SimPad PLUS or LLEAP
- 7) Bilateral - with SimPad Classic, SimPad PLUS or LLEAP
- 8) optional headset/mic - with SimPad Classic, SimPad PLUS or LLEAP
- 9) Anterior and Posterior - with SimPad Classic, SimPad PLUS or LLEAP
- 10) colostomy, ileostomy, cystostomy
- 11) optional SimPad PLUS System
- 12) with interchangeable pupil set
- 13) requires optional Breast Exam Module
- 14) requires optional Mastectomy Module
- 15) requires optional Fundus Skills & Assessment Module
- 16) simulated optional Patient Monitor
- 17) Optional
- 18) Optional - with SimPad Classic, SimPad PLUS or LLEAP
- 19) requires optional SimPad system
- 20) antecubital fossa & dorsal; median, basilic & cephalic
- 21) Requires optional IV arm
- 22) light, tan and brown skin tones
- 23) Medium