



Website: www.dynarex.com
Website: www.thcnet.com

**Fax:** (845) 365-8201

Reorder #'s 1140, 1141, 1143, 1145, 1147

Reviewed on 4/4/13

## **MATERIAL SAFETY DATA SHEET**

# 1. Product and Company Identification

Product/Trade Name: Petrolatum Ointment U.S.P

Manufacturer's Name: Dynarex Corporation

Manufacturer's Address: 10 Glenshaw Street

Orangeburg, NY 10962

Information & Emerg. Phone No.: (888)-DYNAREX (Mon - Fri)

At other times, contact the local Poison Control Center

Chemical Name: Mixture

# 2. Composition/Information on Ingredients

Component NameCAS#/CodesConcentrationRisk PhrasesSafety PhrasesWhite Petrolatum8009-03-08100%R11, R36S2, S16, S25,S26

The above ingredients are raw materials that undergo processing and filling into a three layer sealed laminate. Upon completion of the manufacturing process the material does not possess the occupational health risks associated with the exposure to any of the ingredients above and the finished product.

# 3. Hazards Identification

#### **NFPA Ratings**

NFPA Code for Flammability - 1

NFPA Code for Health - 0

NFPA Code for Reactivity - 0

NFPA Code for Special Hazards - None

## **HMIS Ratings**

HMIS Code for Flammability - 1

HMIS Code for Health - 0

HMIS Code for Reactivity - 0

HMIS Code for Personal Protection - See Section 8





**Fax:** (845) 365-8201

Website: www.dynarex.com
Website: www.thcnet.com

#### **EU Main Hazards**

Flammable.

Irritating to eyes.

## **Routes of Entry**

- Absorption - Eye contact - Ingestion - Inhalation - Skin contact

#### **Carcinogenic Status**

Not considered carcinogenic by NTP, IARC, and OSHA.

# **Target Organs**

- Skin - Eye - Respiratory System

# **Health Effects - Eyes**

The ointment will cause conjunctiva irritation and possibly corneal damage.

#### Health Effects - Skin

No effect on skin, for prolong and short term use. The ointment is for external use only.

May rarely cause allergic reaction in sensitive individuals.

## **Health Effects - Ingestion**

Swallowing may have the following effects:

- irritation of mouth, throat and digestive tract

A large dose may have the following effects:

- Irritation of the stomach and vomiting. May produce mild laxative effect.

#### **Health Effects - Inhalation**

Exposure to vapor may have the following effects:

- irritation of nose, throat and respiratory tract

Exposure to vapor at high concentrations may have the following effects:

- headache





Website: www.dynarex.com
Website: www.thcnet.com

**Fax:** (845) 365-8201

## 4. First Aid Measures

#### Eyes:

Immediately flood the eye with plenty of water for at least 15 minutes, holding the eye open. Obtain medical attention if soreness or redness persists.

## Skin

Immediately flood the skin with large quantities of water, preferably under a shower. Remove contaminated clothing and continue washing. Contaminated clothing should be washed or drycleaned before re-use. Obtain medical attention if blistering occurs or redness persists.

#### Ingestion

Give victim drink 1-3 glasses of water to dilute stomach contents Induce vomiting.

#### Inhalation

Remove from exposure to fresh air. If there is difficulty in breathing, give oxygen. Obtain medical attention immediately.

# 5. Fire Fighting Measures

#### **Extinguishing Media**

Use foam, dry chemical, sand or carbon dioxide; do not use water jet directly. Keep containers and surroundings cool with water spray.

# **Unusual Fire and Explosion Hazards**

Cool endangered containers externally with water.

#### Hazard from combustion products

Nitrogen oxides, Carbon Monoxide, Carbon Dioxide, Soot Organic substances.

# **Protective Equipment for Fire-Fighting**

Do not inhale combustion gases. Wear full protective clothing and self-contained breathing apparatus.





Website: www.dynarex.com
Website: www.thcnet.com

**Fax:** (845) 365-8201

# 6. Accidental Release Measures

Clean area with (hot) water and detergent. Note that spillage will become slippery Avoid contact with hot products Contain and absorb using earth, sand or other inert material. Transfer into suitable containers for recovery or disposal. Wear appropriate protective clothing.

**Environment Precaution:** Do not allow spilled product to enter drains, surface water or ground water. Do not allow product to enter soil/ subsurface. Large spillages must be notified to appropriate authority,

# 7. Handling and Storage

Special precautions are not necessary if use correctly. Keep in tinted close light protecting containers

Use in well ventilated area. Use local exhaust ventilation. Avoid inhaling vapour. Avoid contact with eyes and clothing. Store away from sources of heat or ignition. Storage area should be: - cool - dry - well ventilated - away from incompatible materials (strong oxidizing agents) and heat.

#### 8. Exposure Controls/Personal Protection

## **Occupational Exposure Standards**

There are no specific occupational exposure limits known for product or its ingredients.

# **Engineering Control Measures**

Good general room ventilation is expected to be adequate to control airborne levels.

#### **Respiratory Protection**

Not required under normal conditions and appropriate ventilation.

#### **Hand Protection**

Wear oil resistant gloves made from PVC or nitrile.

#### **Eye Protection**

Chemical goggles or safety glasses with side shields

#### **Body Protection**

Wear oil resistant protecting clothing, overall or apron





Telephone: (845) 365-8200 \* (888) DYNAREX \* Toll Free:

Fax: (845) 365-8201 Website: www.dynarex.com Website: www.thcnet.com

# **Industrial Hygiene**

Do not eat, drink, smoke or sniff during work.

# 9. Physical and Chemical Properties

**Physical State Ointment** 

Color Pale yellow to White Odor Light Lavender Oil odor

pН No data. **Boiling Range/Point** 230°C Melting Point (°C/F) 45°C -60°C Flash Point (PMCC) No data

**Density** 0.780-0.880 gram/ml

Practically insoluble in water. Solubility in Water

No data.

Vapor Density (Air = 1) No data Viscosity (cSt) No data

# 10. Stability and Reactivity

## **Stability**

Stable under normal conditions.

#### **Conditions to Avoid**

**Vapor Pressure** 

- Heat - High temperatures - Static discharge - Exposure to direct sunlight - Mechanical shock

#### **Materials to Avoid**

- Strong oxidizing agents

## **Hazardous Polymerization**

Will not occur.

## **Hazardous Decomposition Products**

- oxides of carbon





200 \* **Website:** www.dynarex.com REX \* **Website:** www.thcnet.com

**Fax:** (845) 365-8201

## 11. Toxicological Information

## **Acute Toxicity**

**Skin:** Rarely cause allergic reaction in sensitive individuals. **Eyes:** Rarely cause allergic reaction in sensitive individuals.

Inhalation: No hazards under normal conditions of use.

# **Chronic Toxicity/Carcinogenicity**

This product is not expected to cause long term adverse health effects.

This product is non carcinogenic.

When swallowed: there is no damage to health.

# Genotoxicity

This product is not expected to cause any mutagenic effects.

# Reproductive/Developmental Toxicity

This product is not expected to cause reproductive or developmental health effects.

## 12. Ecological Information

## **Mobility**

Product is insoluble in water. It can be removed from water mechanically in purification plant.

## Persistence/Degradability

According to the test criteria of OECD the product is not readily biodegradable but inherently biodegradable.

#### **Bio-accumulation**

The product has no bio accumulation.

#### **Ecotoxicity**

**Fish toxicity:** Product is not fish toxic.





\* Website: www.dynarex.com( \* Website: www.thcnet.com

**Fax:** (845) 365-8201

## 13. <u>Disposal Considerations</u>

Dispose of in accordance with all applicable local and national regulations. Labels should not be removed from containers until they have been cleaned. Do not cut, puncture or weld on or near to the container. Do not incinerate closed containers. Empty containers may contain hazardous residues. Dispose of containers with care.

## 14. Transport Information

**DOT CFR 172.101 Data** Not Regulated per 49 CFR 173.4 Small Quantity Exemption

Proper Shipping Name None

Class None
UN Number None
Packing Group None
Classification for AIR None

Transport(IATA):

## 15. Regulatory Information

#### **EU Label Information**

Classification and labeling have been performed according to EU directives 67/548/EEC and 99/45/EC including amendments.

#### **EU Hazard Symbol and Indication of Danger**

F - Highly flammable; Xi - Irritant

# R phrases

R11 - Highly flammable.

R36 - Irritating to eyes.

## S phrases

S2 - Keep out of reach of children.

S16 - Keep away from sources of ignition - No smoking.

S25 - Avoid contact with eyes.

S26 - In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

# US REGULATIONS(Federal, State) and INTERNATIONAL CHEMICAL REGISTRATION LAWS TSCA Listing





Telephone: (845) 365-8200 \*
Toll Free: (888) DYNAREX \*
Fax: (845) 365-8201

\* Website: www.dynarex.com\* Website: www.thcnet.com

This product contains ingredients that are listed on or exempt from listing on the EPA Toxic Substance Control Act Chemical Inventory.

# **EINECS Listing**

All ingredients in this product are listed on the European Inventory of Existing Commercial Chemical Substances (EINECS) or are exempt from listing.

#### **DSL** (Canadian) Listing

All ingredients in this product are listed on the Domestic Substance List (DSL).

## MA Right To Know Law

All components have been checked for inclusion on the Massachusetts Substance List (MSL).

Those components present at or above the de minimal concentration include

# **PA Right To Know Law**

This product contains the following chemicals found on the Pennsylvania Hazardous Substance List

## NJ Right To Know Law

This product contains the following chemicals found on the NJ Right To Know Hazardous Substance List:

## **California Proposition 65**

This product does not contain materials which the State of California has found to cause cancer, birth defects or other reproductive harm.

#### SARA Title III Sect. 302 (EHS)

This product does not contain any chemicals subject to SARA Title III Section 302.

#### SARA Title III Sect. 304

This product does not contain any chemicals subject to SARA Title III Section 304.

# SARA Title III Sect. 311/312 Categorization

Immediate (Acute) Health Hazard

Delayed (Chronic) Health Hazard

Flammable





Fax: (845) 365-8201 Website: www.dynarex.com Website: www.thcnet.com

#### SARA Title III Sect. 313

This product does not contain a chemical which is listed in Section 313 at or above the minimum concentrations.

#### **Abbreviations**

N/A: Denotes no applicable information found or available

CAS#: Chemical Abstracts Service Number

ACGIH: American Conference of Governmental Industrial Hygienists

OSHA: Occupational Safety and Health Administration

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit STEL: Short Term Exposure Limit NTP: National Toxicology Program

IARC: International Agency for Research on Cancer

EU: European Union

R: Risk S: Safety

#### Disclaimer:

The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and Dynarex Corp. assumes no legal responsibility or liability resulting from its use.

This MSDS has been compiled in accordance with - EC Directive 91/155/EC - OSHA's Hazcom Standard (29 CFR 1910.1200).

Use product for intended purpose only, within its expiration date.