

Material Safety Data Sheet

MADACIDE-1

DISTRIBUTED BY:

MADA MEDICAL PRODUCTS, INC.

625 WASHINGTON AVENUE,
CARLSTADT, NJ 07072**DATE PREPARED:** 09/23/11**REVISION NUMBER:** NEW**EMERGENCY #:** 201-460-0454**PRODUCT NUMBER:** 7008, 7009, 7010, 7011**SUPPLIERS PHONE #:** 201-460-0454

SECTION I - IDENTIFICATION

PRODUCT NAME:MADACIDE-1**APPLICATION:**RTU DISINFECTANT

SECTION II - HAZARDOUS INGREDIENTS

HAZARDOUS INGREDIENT	CAS NUMBER	PEL
n-Alkyl (C14, 60%; C16, 30%; C12, 5%; C18, 5%) dimethyl benzyl ammonium chloride	68391-01-5	
n-Alkyl (C12, 68%; C14, 32%) dimethyl ethylbenzyl ammonium chloride	68956-79-6	
Dipropylene glycol monomethyl ether	34590-94-8	

SECTION III - PHYSICAL DATA

APPEARANCE: CLEAR COLORLESS LIQUID/MINTY ODOR**BOILING POINT:** 212 DEGREES F**VAPOR PRESSURE:** NA**VAPOR DENSITY (AIR=1):** NA (Air = 1)**EVAPORATION RATE:** NA**SPECIFIC GRAVITY:** 1.0036**FLAMMABLE LIMITS:** NONCOMBUSTIBLE**SOLUBILITY IN WATER:** COMPLETE**VOLATILITY INCLUDING WATER:** NA

SECTION IV - FIRE AND EXPLOSION DATA

FLASHPOINT: >200 F (PMCC)**EXTINGUISHING MEDIA:** NOT FLAMMABLE**UPPER FLAMMABLE LIMIT (% BY VOL):** NA**LOWER FLAMMABLE LIMIT (% BY VOL):** NA**AUTO-IGNITION TEMPERATURE:** NONE KNOWN**SPECIAL FIRE FIGHTING PROCEDURES:** NONE**UNUSUAL FIRE AND EXPLOSION HAZARDS:** UNKNOWN

Material Safety Data Sheet

MADACIDE-1

SECTION V - REACTIVITY DATA

STABILITY CONDITIONS TO AVOID: NONE
HAZARDOUS POLYMERIZATION:..... WILL NOT OCCUR
HAZARDOUS DECOMPOSITION PRODUCTS: OXIDES OF CARBON. OXIDES OF NITROGEN.
INCOMPATIBILITY: ANIONIC SURFACTANTS, STRONG OXIDIZING
AGENTS
REACTIVITY IN WATER: NONE
CONDITIONS TO AVOID: NONE

SECTION VI - HEALTH DATA

THRESHOLD LIMIT VALUE: NONE
OSHA: NONE
OSHA PEL:..... NONE SET
LISTED CARCINOGEN: NO
MEDICAL CONDITION AGGRAVATED BY EXP: NONE KNOWN
CHRONIC OVER EXPOSURE EFFECTS: REPEATED SKIN CONTACT MAY CAUSE
IRRITATION
IARC MONOGRAPHS: NO
NATIONAL TOXICOLOGY PROGRAM: NO
INHALATION:..... REMOVE TO FRESH AIR. CALL A PHYSICIAN
INGESTION: DO NOT INDUCE VOMITING, DRINK LOTS OF
WATER, CALL A PHYSICIAN
EYES: SPLASH PROOF GOGGLES
SKIN (DERMAL): WASH WITH SOAP & WATER

SECTION VII FIRST AID

BREATHING (INHALATION): REMOVE TO FRESH AIR. CALL PHYSICIAN.
SWALLOWING (INGESTION): DRINK LARGE AMOUNTS OF WATER. DO NOT
INDUCE VOMITING. CALL PHYSICIAN.
EYES: FLUSH WITH WATER FOR 15 MINS. CALL A
PHYSICIAN.
SKIN (DERMAL):..... WASH WITH SOAP & WATER
NOTE TO PHYSICIAN: NONE

SECTION VIII EMPLOYEE PROTECTION

RESPIRATORY PROTECTION:..... NONE NEEDED IN NORMAL USE
LOCAL EXHAUST: ADEQUATE
MECHANICAL EXHAUST: NONE
SPECIAL: NONE

Material Safety Data Sheet

MADACIDE-1

OTHERS: NA
PROTECTIVE CLOTHING: NONE NEEDED IN NORMAL USAGE.
EYE PROTECTION: SPLASH PROOF GOGGLES
PROTECTIVE GLOVES: RUBBER OR PLASTIC

SECTION IX - SPILL AND DISPOSAL DATA

SPILL: SMALL SPILLS FLUSH TO SEWER.

WASTE DISPOSAL: OBSERVE ALL FEDERAL, STATE AND LOCAL REGULATIONS REGARDING DISPOSAL.

OTHER PRECAUTIONS: AVOID PROLONGED SKIN CONTACT

SECTION X - OTHER REGULATORY INFORMATION

HMIS HEALTH: 1
HMIS FLAMMABILITY: 0
HMIS PHYSICAL HAZARD: 0
HMIS OTHER:..... B

SECTION XI - PRECAUTIONARY LABEL STATEMENTS

WARNING: DANGER

SECTION XII - ADDITIONAL INFORMATION

ADDITIONAL: NONE